



PGR Antibody, FITC conjugated

Product Code	CSB-PA11385C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06401
Immunogen	Recombinant Human Progesterone receptor protein (4-232AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Progesterone receptor (PR) (Nuclear receptor subfamily 3 group C member 3), PGR, NR3C3
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PGR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.