



RHBDL2 Antibody, FITC conjugated

Product Code	CSB-PA019672LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NX52
Immunogen	Recombinant Human Rhomboid-related protein 2 protein (12-132AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rhomboid-related protein 2 (RRP2) (EC 3.4.21.105) (Rhomboid-like protein 2) [Cleaved into: Rhomboid-related protein 2, N-terminal fragment (NTF); Rhomboid-related protein 2, C-terminal fragment (CTF)], RHBDL2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RHBDL2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.