




CYP17A1 Antibody

Product Code	CSB-PA006392EA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05093
Immunogen	Recombinant Human Steroid 17-alpha-hydroxylase/17,20 lyase protein (1-508AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Steroid 17-alpha-hydroxylase/17,20 lyase (EC 1.14.14.19) (17-alpha-hydroxyprogesterone aldolase) (EC 1.14.14.32) (CYPXVII) (Cytochrome P450 17A1) (Cytochrome P450-C17) (Cytochrome P450c17) (Steroid 17-alpha-monooxygenase), CYP17A1, CYP17 S17AH
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	CYP17A1

Image



Western blot

All lanes: CYP17A1 antibody at 2µg/ml

Lane 1: HeLa whole cell lysate


Lane 2: A549 whole cell lysate

Secondary


Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 58 kDa


Observed band size: 58 kDa



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA006392EA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA006392EA01HU at dilution of 1:100



Immunofluorescent analysis of A549 cells using CSB-PA006392EA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.