



IFNA14 Antibody, Biotin conjugated

Product Code	CSB-PA011035ED01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01570
Immunogen	Recombinant Human Interferon alpha-14 protein (24-189AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Relevance	Produced by macrophages, IFN-alpha have antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase.
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon alpha-14 (IFN-alpha-14) (Interferon alpha-H) (LeIF H) (Interferon lambda-2-H), IFNA14
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IFNA14