



mela Antibody, Biotin conjugated

Product Code	CSB-PA523525HD01BRG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O34645
Immunogen	Recombinant Bacillus subtilis Alpha-galactosidase protein (1-432AA)
Raised In	Rabbit
Species Reactivity	Bacillus subtilis
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-galactosidase (EC 3.2.1.22) (Melibiase), mela
Immunogen Species	Bacillus subtilis
Research Area	Others
Target Names	mela
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.