




ANGPTL8 Antibody

Product Code	CSB-PA757793LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6UXH0
Immunogen	Recombinant Human Angiopoietin-like protein 8 protein (22-198AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Angiopoietin-like protein 8 (Betatrophin) (Lipasin) (Refeeding-induced fat and liver protein), ANGPTL8, C19orf80 RIFL
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ANGPTL8

Image



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA757793LA01HU at dilution of 1:100