



fim2 Antibody, Biotin conjugated

Product Code	CSB-PA356626LD01BUA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05788
Immunogen	Recombinant Bordetella pertussis Serotype 2 fimbrial subunit protein (27-207AA)
Raised In	Rabbit
Species Reactivity	Bordetella pertussis
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serotype 2 fimbrial subunit, fim2
Immunogen Species	Bordetella pertussis (strain Tohama I / ATCC BAA-589 / NCTC 13251)
Research Area	Others
Target Names	fim2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.