



fim2 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA356626LA01BUA |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P05788 |
| Immunogen | Recombinant Bordetella pertussis Serotype 2 fimbrial subunit protein (27-207AA) |
| Raised In | Rabbit |
| Species Reactivity | Bordetella pertussis |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Serotype 2 fimbrial subunit, fim2 |
| Immunogen Species | Bordetella pertussis (strain Tohama I / ATCC BAA-589 / NCTC 13251) |
| Research Area | Others |
| Target Names | fim2 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |