



# Cd74 Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA004956LD01MO
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P04441
<b>Immunogen</b>	Recombinant Mouse H-2 class II histocompatibility antigen gamma chain protein (56-279AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	H-2 class II histocompatibility antigen gamma chain (Ia antigen-associated invariant chain) (Ii) (MHC class II-associated invariant chain) (CD antigen CD74), Cd74, Ii
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Research Area</b>	Immunology
<b>Target Names</b>	Cd74