





Adck2 Antibody

Product Code	CSB-PA761222LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6NSR3
Immunogen	Recombinant Mouse Uncharacterized aarF domain-containing protein kinase 2 protein (200-609AA)
Raised In	Rabbit
Species Reactivity	Mouse, Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Uncharacterized aarF domain-containing protein kinase 2 (EC 2.7.11.-), Adck2
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Adck2


Image



Western blot
 All lanes: Adck2 antibody at 2µg/ml + K562 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/15000 dilution
 Predicted band size: 69 kDa
 Observed band size: 69 kDa



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA761222LA01MO at dilution of 1:100



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA761222LA01MO at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.