





# Adck2 Antibody

<b>Product Code</b>	CSB-PA761222LA01MO
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q6NSR3
<b>Immunogen</b>	Recombinant Mouse Uncharacterized aarF domain-containing protein kinase 2 protein (200-609AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Mouse, Human
<b>Tested Applications</b>	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Uncharacterized aarF domain-containing protein kinase 2 (EC 2.7.11.-), Adck2
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Research Area</b>	Others
<b>Target Names</b>	Adck2


## Image



Western blot  
All lanes: Adck2 antibody at 2µg/ml + K562 whole cell lysate  
Secondary  
Goat polyclonal to rabbit IgG at 1/15000 dilution  
Predicted band size: 69 kDa  
Observed band size: 69 kDa



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA761222LA01MO at dilution of 1:100



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA761222LA01MO at dilution of 1:100