



Hemoglobin Antibody

Product Code	CSB-PA00500E0Gp
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Rabbit Hemoglobin (Native Protein)
Raised In	Guinea pig
Species Reactivity	Rabbit
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Haemoglobin
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.