



Mup2 Antibody

Product Code	CSB-PA319316YA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11589
Immunogen	Recombinant Mouse Major urinary protein 2 protein (19-180AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major urinary protein 2 (MUP 2), Mup2
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Mup2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.