



yoaJ Antibody, Biotin conjugated

Product Code	CSB-PA522026YD01BRJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O34918
Immunogen	Recombinant Bacillus subtilis Expansin-YoaJ protein (26-232AA)
Raised In	Rabbit
Species Reactivity	Bacillus subtilis
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Expansin-YoaJ (EXLX1), yoaJ
Immunogen Species	Bacillus subtilis (strain 168)
Research Area	Others
Target Names	yoaJ
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.