



Mup11 Antibody

Product Code	CSB-PA361419YA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04938
Immunogen	Recombinant Mouse Major urinary proteins 11 and 8 protein (1-151AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major urinary protein 11, Mup11, Mup9
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Mup11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.