



# Spp1 Antibody, Biotin conjugated

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | CSB-PA022603YD01MO  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.   |
| <b>Uniprot No.</b>         | P10923  |
| <b>Immunogen</b>           | Recombinant Mouse Osteopontin protein (17-290AA)  |
| <b>Raised In</b>           | Rabbit  |
| <b>Species Reactivity</b>  | Mouse   |
| <b>Tested Applications</b> | ELISA   |
| <b>Form</b>                | Liquid  |
| <b>Conjugate</b>           | Biotin  |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4  |
| <b>Purification Method</b> | >95%, Protein G purified  |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Alias</b>               | Osteopontin (2AR) (Bone sialoprotein 1) (Calcium oxalate crystal growth inhibitor protein) (Early T-lymphocyte activation 1 protein) (Minopontin) (Secreted phosphoprotein 1) (SPP-1), Spp1, Eta-1 Op Spp-1 |
| <b>Immunogen Species</b>   | Mus musculus (Mouse)  |
| <b>Research Area</b>       | Cardiovascular  |
| <b>Target Names</b>        | Spp1  |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.   |