



Spp1 Antibody, Biotin conjugated

Product Code	CSB-PA022603YD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10923
Immunogen	Recombinant Mouse Osteopontin protein (17-290AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Osteopontin (2AR) (Bone sialoprotein 1) (Calcium oxalate crystal growth inhibitor protein) (Early T-lymphocyte activation 1 protein) (Minopontin) (Secreted phosphoprotein 1) (SPP-1), Spp1, Eta-1 Op Spp-1
Immunogen Species	Mus musculus (Mouse)
Research Area	Cardiovascular
Target Names	Spp1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.