



# CSF1R Antibody

<b>Product Code</b>	CSB-PA006044LA01CA
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P13369
<b>Immunogen</b>	Recombinant Felis silvestris catus Macrophage colony-stimulating factor 1 receptor protein (20-514AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Felis silvestris catus
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Macrophage colony-stimulating factor 1 receptor (CSF-1 receptor) (CSF-1-R) (CSF-1R) (M-CSF-R) (EC 2.7.10.1) (Proto-oncogene c-Fms) (CD antigen CD115), CSF1R, FMS
<b>Immunogen Species</b>	Felis catus (Cat) (Felis silvestris catus)
<b>Research Area</b>	Others
<b>Target Names</b>	CSF1R
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.