



Ms4a1 Antibody, Biotin conjugated

Product Code	CSB-PA015007LD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19437
Immunogen	Recombinant Mouse B-lymphocyte antigen CD20 protein (111-291AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	B-lymphocyte antigen CD20 (B-cell differentiation antigen Ly-44) (Lymphocyte antigen 44) (Membrane-spanning 4-domains subfamily A member 1) (CD antigen CD20), Ms4a1, Cd20 Ly-44 Ms4a2
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	Ms4a1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.