



TLR2 Antibody, HRP conjugated

Product Code	CSB-PA023601LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60603
Immunogen	Recombinant Human Toll-like receptor 2 protein (19-588AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Toll-like receptor 2 (Toll/interleukin-1 receptor-like protein 4) (CD antigen CD282), TLR2, TIL4
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TLR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.