



TLR2 Antibody

Product Code	CSB-PA023601LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60603
Immunogen	Recombinant Human Toll-like receptor 2 protein (19-588AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Toll-like receptor 2 (Toll/interleukin-1 receptor-like protein 4) (CD antigen CD282), TLR2, TIL4
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TLR2