



# TLR2 Antibody

<b>Product Code</b>	CSB-PA023601LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O60603
<b>Immunogen</b>	Recombinant Human Toll-like receptor 2 protein (19-588AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Toll-like receptor 2 (Toll/interleukin-1 receptor-like protein 4) (CD antigen CD282), TLR2, TIL4
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Immunology
<b>Target Names</b>	TLR2
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.