



gG Antibody, Biotin conjugated

Product Code	CSB-PA361844HD01HWY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06484
Immunogen	Recombinant Human herpesvirus 1 Envelope glycoprotein G protein (25-189AA)
Raised In	Rabbit
Species Reactivity	Human herpesvirus 1
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Envelope glycoprotein G (gG) (gG-1), gG
Immunogen Species	Human herpesvirus 1 (strain 17) (HHV-1) (Human herpes simplex virus 1)
Research Area	Others
Target Names	gG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.