



ADIPOQ Antibody

Product Code	CSB-PA661101LA01BO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3Y5Z3
Immunogen	Recombinant Bovine Adiponectin protein (18-240AA)
Raised In	Rabbit
Species Reactivity	Bovine
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Adiponectin (30 kDa adipocyte complement-related protein) (Adipocyte complement-related 30 kDa protein) (ACRP30) (Adipocyte, C1q and collagen domain-containing protein) (Adipose most abundant gene transcript 1 protein) (apM-1), ADIPOQ, ACRP30 APM1
Immunogen Species	Bos taurus (Bovine)
Research Area	Others
Target Names	ADIPOQ
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.