



Flii Antibody

Product Code	CSB-PA884922LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9JJ28
Immunogen	Recombinant Mouse Protein flightless-1 homolog protein (495-827AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein flightless-1 homolog, Flii, Fli1 Fliih
Immunogen Species	Mus musculus (Mouse)
Research Area	Signal Transduction
Target Names	Flii
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.