



GT-1 Antibody

Product Code	CSB-PA861892LA01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9FX53
Immunogen	Recombinant Arabidopsis thaliana Trihelix transcription factor GT-1 protein (1-406AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Trihelix transcription factor GT-1 (Trihelix DNA-binding protein GT-1), GT-1
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Research Area	Others
Target Names	GT-1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.