




AHCY Antibody

Product Code	CSB-PA001474LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23526
Immunogen	Recombinant Human Adenosylhomocysteinase protein (2-432AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Adenosylhomocysteinase (AdoHcyase) (EC 3.3.1.1) (S-adenosyl-L-homocysteine hydrolase), AHCY, SAHH
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	AHCY

Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA001474LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.