



CETN2 Antibody

Product Code	CSB-PA005265HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41208
Immunogen	Recombinant Human Centrin-2 protein (1-172AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Centrin-2 (Caltractin isoform 1), CETN2, CALT CEN2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CETN2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.