





FBXL18 Antibody

Product Code	CSB-PA008469LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96ME1
Immunogen	Recombinant Human F-box/LRR-repeat protein 18 protein (1-365AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	F-box/LRR-repeat protein 18 (F-box and leucine-rich repeat protein 18), FBXL18, FBL18
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FBXL18


Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA008469LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human lymphoid tissue using CSB-PA008469LA01HU at dilution of 1:100



Immunofluorescent analysis of A549 cells using CSB-PA008469LA01HU at dilution of 1:100 and Cy3-conjugated Goat Anti-Rabbit IgG

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.