



UCP1 Antibody, HRP conjugated

Product Code	CSB-PA025554LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25874
Immunogen	Recombinant Human Mitochondrial brown fat uncoupling protein 1 protein (2-307AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitochondrial brown fat uncoupling protein 1 (UCP 1) (Solute carrier family 25 member 7) (Thermogenin), UCP1, SLC25A7 UCP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	UCP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.