



FOXP1 Antibody, FITC conjugated

Product Code	CSB-PA008841EC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H334
Immunogen	Recombinant Human Forkhead box protein P1 protein (1-114AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Relevance	Transcriptional repressor. It plays an important role in the specification and differentiation of lung epithelium. Can act with CTBP1 to synergistically repress transcription but CTBP1 is not essential. Essential transcriptional regulator of B-cell development
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Forkhead box protein P1 (Mac-1-regulated forkhead) (MFH), FOXP1
Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	FOXP1