



Amh Antibody, Biotin conjugated

Product Code	CSB-PA001666YD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P27106
Immunogen	Recombinant Mouse Muellierian-inhibiting factor protein (450-552AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Muellerian-inhibiting factor (Anti-Muellerian hormone) (AMH) (Muellerian-inhibiting substance) (MIS), Amh
Immunogen Species	Mus musculus (Mouse)
Research Area	Signal Transduction
Target Names	Amh