



Mybpc3 Antibody

Product Code	CSB-PA015267YA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56741
Immunogen	Recombinant Rat Myosin-binding protein C, cardiac-type protein (645-864AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myosin-binding protein C, cardiac-type (Cardiac MyBP-C) (C-protein, cardiac muscle isoform), Mybpc3
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Cardiovascular
Target Names	Mybpc3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.