



RPL35A Antibody, HRP conjugated

Product Code	CSB-PA020248LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18077
Immunogen	Recombinant Human 60S ribosomal protein L35a (1-110AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	60S ribosomal protein L35a (Cell growth-inhibiting gene 33 protein) (Large ribosomal subunit protein eL33), RPL35A
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RPL35A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.