





RPL35A Antibody

Product Code	CSB-PA020248LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18077
Immunogen	Recombinant Human 60S ribosomal protein L35a protein (1-110AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	60S ribosomal protein L35a (Cell growth-inhibiting gene 33 protein) (Large ribosomal subunit protein eL33), RPL35A
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RPL35A

Image



Immunohistochemistry of paraffin-embedded human lymphoid tissue using CSB-PA020248LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human spleen tissue using CSB-PA020248LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.