



PODXL Antibody, HRP conjugated

Product Code	CSB-PA13447B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00592
Immunogen	Recombinant Human Podocalyxin protein (32-458AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Podocalyxin (GCTM-2 antigen) (Gp200) (Podocalyxin-like protein 1) (PC) (PCLP-1), PODXL, PCLP PCLP1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PODXL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.