



NUP210 Antibody, FITC conjugated

Product Code	CSB-PA016195LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TEM1
Immunogen	Recombinant Human Nuclear pore membrane glycoprotein 210 protein (28-238AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear pore membrane glycoprotein 210 (Nuclear pore protein gp210) (Nuclear envelope pore membrane protein POM 210) (POM210) (Nucleoporin Nup210) (Pore membrane protein of 210 kDa), NUP210, KIAA0906
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	NUP210
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.