



BMP2 Antibody, FITC conjugated

Product Code	CSB-PA09419C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12643
Immunogen	Recombinant Human Bone morphogenetic protein 2 protein (283-396AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Bone morphogenetic protein 2 (BMP-2) (Bone morphogenetic protein 2A) (BMP-2A), BMP2, BMP2A
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	BMP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.