



Atp5f1b Antibody

Product Code	CSB-PA002350LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56480
Immunogen	Recombinant Mouse ATP synthase subunit beta, mitochondrial protein (47-529AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATP synthase subunit beta, mitochondrial (ATP synthase F1 subunit beta), Atp5f1b, Atp5b
Immunogen Species	Mus musculus (Mouse)
Research Area	Tags & Cell Markers
Target Names	Atp5f1b
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.