





CCL8 Antibody

Product Code	CSB-PA07629A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P80075
Immunogen	Recombinant Human C-C motif chemokine 8 protein (24-99AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-C motif chemokine 8 (HC14) (Monocyte chemoattractant protein 2) (Monocyte chemotactic protein 2) (MCP-2) (Small-inducible cytokine A8) [Cleaved into: MCP-2(6-76)], CCL8, MCP2 SCYA10 SCYA8
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CCL8

Image



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA07629A0Rb at dilution of 1:100



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA07629A0Rb at dilution of 1:100