



Disintegrin EC3B Antibody

Product Code	CSB-PA305523HA01EAH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P81631
Immunogen	Recombinant Echis carinatus Disintegrin EC3B protein (1-67AA)
Raised In	Rabbit
Species Reactivity	Echis carinatus
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Disintegrin EC3B
Immunogen Species	Echis carinatus
Research Area	Others