



TXNL1 Antibody, Biotin conjugated

Product Code	CSB-PA025384HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43396
Immunogen	Recombinant Human Thioredoxin-like protein 1 protein (2-289AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thioredoxin-like protein 1 (32 kDa thioredoxin-related protein), TXNL1, TRP32 TXL TXNL
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TXNL1