



ROPN1B Antibody, Biotin conjugated

Product Code	CSB-PA020065ED01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BZX4
Immunogen	Recombinant Human Ropporin-1B protein (1-120AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ropporin-1B (Rhophilin-associated protein 1B), ROPN1B
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ROPN1B