



IFNA10 Antibody, Biotin conjugated

Product Code	CSB-PA011033LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01566
Immunogen	Recombinant huamn Interferon alpha-10 protein (24-189AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon alpha-10 (IFN-alpha-10) (Interferon alpha-6L) (Interferon alpha-C) (LeIF C), IFNA10
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IFNA10