





IFNA10 Antibody

Product Code	CSB-PA011033LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01566
Immunogen	Recombinant huamn Interferon alpha-10 protein (24-189AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon alpha-10 (IFN-alpha-10) (Interferon alpha-6L) (Interferon alpha-C) (LeIF C), IFNA10
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IFNA10

Image



Western blot
 All lanes: IFNA10 antibody at 2µg/ml + Mouse kidney tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 22 kDa
 Observed band size: 22 kDa



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA011033LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.