



KRT10 Antibody, HRP conjugated

Product Code	CSB-PA012504LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13645
Immunogen	Recombinant Human Keratin, type I cytoskeletal 10 protein (1-584AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Keratin, type I cytoskeletal 10 (Cytokeratin-10) (CK-10) (Keratin-10) (K10), KRT10, KPP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KRT10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.