



# manA Antibody, Biotin conjugated

Product Code	CSB-PA014754HD01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00946
Immunogen	Recombinant Escherichia coli Mannose-6-phosphate isomerase protein (1-391AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mannose-6-phosphate isomerase (EC 5.3.1.8) (Phosphohexomutase) (Phosphomannose isomerase) (PMI), manA, pmi
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	manA