



manA Antibody, HRP conjugated

Product Code	CSB-PA014754HB01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00946
Immunogen	Recombinant Escherichia coli Mannose-6-phosphate isomerase protein (1-391AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mannose-6-phosphate isomerase (EC 5.3.1.8) (Phosphohexomutase) (Phosphomannose isomerase) (PMI), manA, pmi
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	manA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.