




# manA Antibody

|                            |                                                                                                            |
|----------------------------|------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA014754HA01ENV                                                                                        |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                              |
| <b>Uniprot No.</b>         | P00946                                                                                                     |
| <b>Immunogen</b>           | Recombinant Escherichia coli Mannose-6-phosphate isomerase protein (1-391AA)                               |
| <b>Raised In</b>           | Rabbit                                                                                                     |
| <b>Species Reactivity</b>  | Escherichia coli                                                                                           |
| <b>Tested Applications</b> | ELISA, WB; Recommended dilution: WB:1:1000-1:5000                                                          |
| <b>Form</b>                | Liquid                                                                                                     |
| <b>Conjugate</b>           | Non-conjugated                                                                                             |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4                           |
| <b>Purification Method</b> | >95%, Protein G purified                                                                                   |
| <b>Isotype</b>             | IgG                                                                                                        |
| <b>Clonality</b>           | Polyclonal                                                                                                 |
| <b>Alias</b>               | Mannose-6-phosphate isomerase (EC 5.3.1.8) (Phosphohexomutase) (Phosphomannose isomerase) (PMI), manA, pmi |
| <b>Immunogen Species</b>   | Escherichia coli                                                                                           |
| <b>Research Area</b>       | Others                                                                                                     |
| <b>Target Names</b>        | manA                                                                                                       |

## Image



### Western blot

All lanes: manA antibody at 16μg/ml + DH5α whole cell lysate

### Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 43 kDa

Observed band size: 43 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.