



Defb14 Antibody, Biotin conjugated

Product Code	CSB-PA749254HD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7TNV9
Immunogen	Recombinant Mouse Beta-defensin 14 protein (23-67AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-defensin 14 (BD-14) (mBD-14) (Defensin, beta 14), Defb14
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Defb14
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.